

Certificate of Analysis

Pertuzumab 单抗特异性抗原

Cat. No.: QRE-105
Lot. No.: AG050180602
Expiration Date: May 2021

Storage: The product is shipped with dry ice and should be stored below -20°C upon arrival.

Specificity by ELISA: Conforms.

Protein Purity Assay (SDS-PAGE): no less than 95% pure as determined by SDS-PAGE analysis using Coomassie Blue detection.

Protein Concentration by E280: 1.0 mg/ml

Extinction Coefficient at 0.1% = 0.945

For Research Use Only. Not for diagnostic use.

Suzhou Rhinozyme Biotechnology Co., Ltd.

No. 218, Xinghu Street, Suzhou Industrial Park

Jiangsu Province, China, 215123

For inquiries, contact us at 0086-512-87663137 or techserv@rhinobio.com.

www.rhinobio.com

H. J.
2018年06月01日
QA/QC Manager
Issued on 01 Jun 2018

