

# Certificate of Analysis

## Pertuzumab 单抗特异性抗原

Cat. No.: QRE-105  
Lot. No.: AG050180903  
Expiration Date: Aug. 2021

---

**Storage:** The product is shipped with dry ice and should be stored below -20°C upon arrival.

**Specificity by ELISA:** Conforms.

**Protein Purity Assay (SDS-PAGE) :** no less than 95% pure as determined by SDS-PAGE analysis using Coomassie Blue detection.

**Protein Concentration by E280:** 1.0 mg/ml

**Extinction Coefficient at 0.1% = 0.945**

---

For Research Use Only. Not for diagnostic use.

### Suzhou Rhinozyme Biotechnology Co., Ltd.

No. 218, Xinghu Street ,Suzhou Industrial Park

Jiangsu Province, China,215123

For inquiries, contact us at 0086-512-87663137or techserv@rhinobio.com.

www.rhinobio.com

---

QA/QC Manager  
Issued on 10 Sep 2018



**RHINO BIO**  
瑞安生物科技