

# Certificate of Analysis

## 重组人 ErbB2 抗原

Cat. No.: QRE-103  
Lot. No.: AG030231201  
Expiration Date: Nov 2026

**Storage:** The product is shipped with dry ice and should be stored below -20°C upon arrival.

**Specificity by ELISA:** Conforms.

**Protein Purity Assay (SDS-PAGE) :** no less than 95% pure as determined by SDS-PAGE analysis using Coomassie Blue detection.

**Protein Concentration by E280:** 1.0 mg/ml

**Extinction Coefficient** at 0.1% = 0.953

For Research Use Only. Not for diagnostic use.

### Shanghai RhinoBio Co.,Ltd

99 Ruoshui Road, Suzhou Industrial Park, Jiangsu Province .

For inquiries, contact us at 0086-512-87663137or techserv@rhinobio.com.

www.rhinobio.com

