

# Certificate of Analysis

## Human FcγRI (His 标签)

Cat. No.: QRE-133  
Lot. No.: AG330231102  
Expiration Date: Oct 2024

**Storage:** The product is shipped with dry ice and should be stored below -20°C upon arrival.

**Specificity by ELISA:** Conforms.

**Protein Purity Assay (SDS-PAGE):** no less than 95% pure as determined by SDS-PAGE analysis using Coomassie Blue detection.

**Protein Concentration by E280:** 1.0mg/ml.

**Extinction Coefficient at 0.1% = 1.772.**

For Research Use Only. Not for diagnostic use.

### Shanghai RhinoBio Co.,Ltd

99 Ruoshui Road, Suzhou Industrial Park, Jiangsu Province .  
For inquiries, contact us at 0086-512-87663137or techserv@rhinobio.com.  
www.rhinobio.com

  
QA/QC Manager  
Issued on 17 Nov 2023